powerball® USER MANUAL AUTOSTART™

or visit:

POWERBALLS.COM/START

to learn the perfect starting & operating technique

3 THINGS YOU NEED TO KNOW

Always spin **ONLY**
with your wrist.
The arm **MUST** be
kept 100% still.

If the ball only spins
for 15 seconds and then
stops, the issue is with your
spin technique and not the
ball.

Please watch this short
video:

powerballs.com/start

It will immediately
show you the correct
technique and how to
reach those really high
spin speeds.

If you still have a question, issue, query or problem with your new Powerball® then please email:

help@powerballs.com

Anything at all - just contact us and we'll respond by return to solve it.

Can't start Powerball®? powerballs.com/start

Have a question or query about your Powerball®? help@powerballs.com

Need to chat with us? Call or Whatsapp: +353 83 422 6312

Become a stockist: sales@powerballs.com

facebook.com/powerball

twitter.com/powerball

General Precautions

- Powerball® is an isometric exercise/rehabilitation device which can produce powerful resistance when operated at high speed. Therefore, if in any doubt about your own physical ability, either operate at slow spin speed or consult with your physician before operating this product.
- Powerball® is suitable for all ages. However, because of its potential to generate high resistance, it is recommended that children UNDER 12 are supervised by an adult while using this product to help avoid strain on developing upper limbs.
- Avoid touching the rotor while it is spinning as doing so may cause friction burn and/or pinch the skin.
- Please keep the inner part of your Powerball® clean and avoid any form of lubricant or oil contamination. If accidentally oiled, or if foreign material is ingested into the sphere, your Powerball® may be easily opened and cleaned - please visit powerballs.com for details.
- Always aim to grip Powerball® FIRMLY in the hand each time you spin - doing so will help to extend the lifespan of your Powerball®.
- Try to avoid dropping Powerball® at any time.

For starting, operating & maintenance tips, visit:
POWERBALLS.COM/START

1 WIND...

Use one thumb to hold, one to wind forward...

2 WIND 4 FULL TURNS & STOP

...you'll be able to feel the 'stop' point as the rotor won't wind any further

3 HOLD...

Hold the fully wound rotor in place...

4 RELEASE & SPIN!

...grip the ball firmly in your preferred hand & release the rotor

A

Slow lazy turns of the hand

B

...never spin with the arm; it's ALWAYS with the wrist ONLY!

LCD SPEED METER FUNCTIONS

Total Accumulated Revolutions

Press the 'FCN' button to reveal the 'total revolutions' mode on the display. This mode calculates the accumulated rotor revolutions during a session and records each 100 rotations as one unit - for example 7240 revolutions will display as 72. This is very useful when rehabilitating a damaged upper limb where low speed / long spin sessions are preferable to short, high speed spins.

Real Time Revolutions

Press the 'FCN' button to reveal the 'current revolutions' mode on the display. This counter setting allows real time monitoring of the rotor's RPM (revs per minute) speed and shows the actual rotor speed at any specific time during the spin session.

*** Note: The LCD display is induction powered & only operates once the rotor speed exceeds 1,800rpm**

The Highest Achieved Spin Speed

The 'highest speed' function stores and displays the historical highest spin speed achieved on your Powerball® over the past 48hr period (after which it will be automatically deleted). Press and hold the 'FCN/CLR' button to clear the high speed score manually. This can be very useful for keeping track of strength gains while rehabilitating an injured limb or simply for enjoyable competition among friends.

The Muscular Endurance Mode

Press the 'FCN' button to reveal the 'strength index' mode on the display. This setting will begin to record the number of actual revolutions made during the 30, 60 or 90 second interval immediately after the rotor is put in motion (press the FCN/CLR button to cycle between all three). The display begins counting down and continues recording until the last 5 seconds - whereby the total accumulated revs will appear as it counts down to '0'. This is a superb test of muscular endurance as you'll need to keep the rotor spinning at the highest possible spin speed for a period of 30, 60 or 90 seconds in order to achieve a maximum score.

«« 'Pro' models feature induction powered **multi function digital counters** to precisely measure spin speed performance data. Includes automatic backlit display. They power on once the rotor reaches 1,800rpm

«« Fusion & Blaze Powerball® models feature **magnetic induction powered** colour LEDs built into the rotor which glow beautifully as the ball spins

«« **Auto start mechanism** for fuss free starts every time - simply wind backwards 3-4 turns in the direction of the arrows inscribed on the rotor, release & start spinning

«« The patented 280 Series Powerball® shell casing is **drop resistant** and easy to open for cleaning. It features reversible inner tracks for **2X** performance lifespan

«« All 280 Powerball® rotors are **certified vibration free to 17,000rpm spin speeds** ensuring powerful resistance and a 100% non-impact spin experience

GETTING FASTER WITH YOUR POWERBALL®?

**Submit your high scores to the
international Powerball® scoreboards
& compare your spin speed to
Powerballers around the world, with
great prizes each month for the
fastest spin!**

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All images are protected by trade marks law and/or copyright.
© RPM Sports Limited 2018 - All rights reserved.
Distributed by: RPM Sports Ltd